



OXYGENCONCENTRATORSTORE

FROM AMERICAN MEDICAL SALES & RENTALS

Home Oxygen Concentrator Comparison

Unit Name	Oxygen Flow Options	Continuous Flow Liters	Weight	Dimensions (WxHxD) inches
Inogen at Home 5L				
	Continuous	1 to 5 LPM	18 lbs (or 8.2 kg)	13 x 16.5 x 7 in
OxyGo OxyHome 5L				
	Continuous	0.5 to 5 LPM	41.8 lbs (or 19 kg)	13.4 x 25.6 x 11.8
Drive DeVilbiss 525DS Quiet				
	Continuous	0.5 to 5 LPM	36 lbs (or 16.3 kg)	13.5 x 24.5 x 12 in
Drive DeVilbiss 1025DS				
	Continuous	2 to 10 LPM	42 lbs (or 19 kg)	13.5 x 24.5 x 12 in
CAIRE Intensity 10				
	Continuous	2 to 10 LPM	58 lbs (or 26.3 kg)	16.5 x 27.5 x 14.5 in



Unit Name	Noise	Humidifier Bottle Compatible	CPAP Machine Compatible	Power Consumption	Max Operating Altitude	Warranty
Inogen at Home 5L						
	Less than 50 dBA	Yes	Yes	Less than 275 W	0 to 10,000 feet	Three-Year Manufacturer's Warranty
OxyGo OxyHome 5L						
	Less than 45 dBA	Yes	Yes	Less than 350 W	0 to 10,000 feet	Three-Year Manufacturer's Warranty
Drive DeVilbiss 525DS Quiet						
	46 dBA	Yes	Yes	Less than 310 W	0 to 13,123 feet	Three-Year Manufacturer's Warranty
Drive DeVilbiss 1025DS						
	Less than 59 dBA	Yes	Yes	639 W	0 to 5,000 feet	Three-Year Manufacturer's Warranty
CAIRE Intensity 10						
	58 dBA	Yes	Yes	600 W	-1,250 to 10,000 feet	Three-Year Manufacturer's Warranty