

P2-S3 PORTABLE OXYGEN CONCENTRATOR



ITEM #	DESCRIPTION
P2-S3	Portable Oxygen Concentrator

- Compact & Lightweight.
- High Efficiency, Low Noise. Excellent Heat Dissipation.
- High Adsorption Sieve Bed.
- Long Lasting Lithium-Ion Battery - Up to 7.2 Hours Battery Life.
- Single Battery Standard. Available Optional Double and Triple Battery.
- Includes AC/DC Car Charger.
- Conforms to All Applicable FAA Acceptance Criteria for POC Carriage and Use on Board Aircraft.

SPECIFICATIONS:

Size	3.4" W x 6.2" D x 6.7" H (with Single Battery)
Weight	3.3 lbs
Warm up time	2 Minutes
Settings	Three Flow Settings (1-3)
Oxygen Concentration	87%-96%
AC Power Supply	100-240v, 50-60hz
DC Power Supply	1-16V DC Output: 19V 6.3A Voltage: 14.4V
Rechargeable Battery	14.4V DC 3.35A MAX Single Battery Rated Capacity: 3,400 mAh Double Battery Rated Capacity: 6,800 mAh Triple Battery Rated Capacity: 10,200 mAh
Display Screen	2"
Operation	Simple Control Functions and Easy-to-Read LCD Display
Use	Designed for 24/7 Use

BATTERY DURATION

Single Battery

Setting 1: 2.5h
Setting 2: 1.7h
Setting 3: 1.2h

Double Battery

Setting 1: 4.8h
Setting 2: 3.5h
Setting 3: 2.5h

Triple Battery

Setting 1: 7.2h
Setting 2: 5h
Setting 3: 3.5h

WARRANTY

Three Years. Sieve Bed: Two Years. Battery: One Year

V1: 2/2/24