

Specifications

Model:	P2
Size:	3.3" W x 8.7" D x 6.3" H
Weight:	4.37 lbs
Warm up time:	2 minutes
Settings:	Five Flow Settings (1-5)
Oxygen Concentration:	87%-96%
Oxygen Flow:	Pulse dose delivery system, 5 flow settings
AC Power Supply:	100-240v, 50-60hz (Auto sensing to allow for worldwide use.)
Rechargeable Battery:	14.4VDC
Battery Rated Capacity:	6.8Ah
Display Screen:	2.8"
Battery:	6800mAh, duration up to 5 hours on setting 1 Approx. 4 hours charge time with AC
Operation:	Simple control functions and easy-to-read LCD display
Use:	Designed for 24/7 use
Warranty:	Three years, excluding battery and sieve bed

Portable Oxygen Concentrator

P2



American Medical Sales and Rentals

15055 E. Hinsdale Drive, Bldg. D Centennial, CO 80112
Phone: +1 (877) 774-9271

www.oxygenconcentratorstore.com
E-mail: sales@amsrco.com



Compact & Lightweight



External Battery Charger*


Conveniently charge a backup battery while the device is in use.

Charging Time: 2.5 hours

* Additional purchase.



User Control



Powerful Compressor



High Efficiency




Low Noise



Excellent Heat
Dissipation



Long Life



Wherever you go Oxygen follows



4.4 lbs
Lightweight



Equivalent
to 5 Liters



High Absorption
Sieve Bed



Rechargeable
Battery
Duration up to 4 hrs



93% O₂
Concentration



Small and
Convenient

