

P2 PORTABLE OXYGEN CONCENTRATOR



ITEM #	DESCRIPTION
P2	1L Portable Oxygen Concentrator

- Compact & lightweight.
- Conveniently charge a backup battery while the device is in use. Charging time 2.5 hours. (additional purchase).
- High efficiency, low noise. Excellent heat dissipation.
- 93% O₂ concentration.
- High absorption sieve bed.
- Equivalent up to 5L.
- Includes AC/DC Car Charger.
- Meet FAA Travel Restrcitions.

FEATURES

Size	8.7" x 3.4" x 6.3"
Weight	4.37lbs
Warm up time	2 minutes
Settings	Five Flow Settings (1-5)
Oxygen Concentrator	87%-96%
AC Power Supply:	100-240v, 50-60hz (auto sensing to allow for worldwide use)
Rechargeable Battery	14.4VDC
Battery Rated Capacity:	6.8Ah
Display Screen	2.8"
Battery	6800mAh, duration up to 5 hours on setting 1. Approx. 4 hours charge time with AC.
Operation	Simple control functions and easy-to-read LCD display
Use	Designed for 24/7 use

WARRANTY

Three years, excluding battery and sieve bed

V3: 4/22/21